USER MANUAL

Bedside Navigator with UCL FISECARE v2.0

Running FISECARE with MotionInput's Bedside Navigation

Step 1: Open Bedside Navigator. The Settings Window will open.

🔏 UCL MotionInput v3.2	Bedside Navigator	×
General Navigation with	i Face Hand	About ? UCL FISECARE
Calibration Default Camera	i OFF i Camera 1 ~	Open FISECARE i OFF
Low Light	i ON	Edit Apps i Open App Settings
Facial Navigation	i Nose O Joystick	Hand Navigation Handedness i OLeft ORight
Method Mouse Speed	i Speech Gestures	
		Launch Cancel

Step 2: Select your preferred settings for navigation.

Select **Hand** if navigating the cursor with either your left or right hand OR select **Face** if navigating the cursor **using your nose**.

Click the **Information buttons** highlighted in **red**, beside each option to learn more about the modes. Including what gestures and speech commands are available to help you navigate FISECARE.

Finally, click the 'Launch' button, shown highlighted in blue to start MotionInput.

The 'Open App Settings' button allows you to edit and personalise apps on the FISECARE interface. If '**Open App Settings'** in **orange** is chosen, then follow Steps 3 onwards.

If the button highlighted in **green** is chosen, then follow Steps 4 to Steps 6.

To Launch UCL FISECARE, make sure '**LAUNCH WITH FISECAR**E' is set to '**YES'**. You may need to install <u>.NET Desktop Runtime 6.0</u>.

🔏 UCL MotionInput v3.2	Bedside Navigator	×
General	i © Face O Hand	About ?
Navigation with		UCL FISECARE
Calibration	i OFF	Open FISECARE i OFF
Default Camera	i Camera 1 ~	Edit Apps i Open App Settings
Low Light	i ON	
Facial Navigation		Hand Navigation
Mode	i Nose OJoystick	Handedness i O Left O Right
Method	i Speech Gestures	
Mouse Speed	i 5	
		Launch Cancel

Step 3: Edit Apps on FISECARE

FISECARE App Settings				×
Example:				
App name: B	edside Navigator	App File Path:	C:\Program Files\Bedside Navigator\Bedside Navigator.	exe
Enter the name	Enter the name of the app		te path of the app starting with C:\	
App 1 Name	Feet Navigation	App 1 File Path	C:\Program Files\Feet Navigation\Feet Navigation.	Select File
App 2 Name	Bedside Navigator	App 2 File Path	C:\Program Files\Bedside Navigator\Bedside Navi	Select File
App 3 Name	None	App 3 File Path	None	Select File
App 4 Name	None	App 4 File Path	None	Select File
App 5 Name	None	App 5 File Path	None	Select File
			Save	Cancel

A new window like the image above will pop up. Click 'Select File' and browse through your files to find the application's executable file (.exe). Then give each application a name.

For example,

FISECARE App Settings						>	< <u> </u>	, ○ Search Internet
Example:								
App name: B	edside Navigator		App File Path: C:\P	Program Files\Bed	side Navigator\Bedside Navigat	or.exe		
Enter the name	of the app	s	elect the absolute pa	ath of the app star	ting with C:\			
App 1 Name	None	A	pp 1 File Path	lone		Select File		
			🦺 Open					×
App 2 Name	None	A	← → ~ ↑ .	« Windows (C:)	Program Files ⇒ Internet Explore	er ~ Č		
			Organise - Ne	w folder			8	- 🔳 🕐
App 3 Name	None	A	📙 Settings	^ Name		Date modified	Туре	Size
			💻 Desktop	📙 en-US		07/12/2019 2:44 PM	File folder	
			OneDrive	📙 image	es	07/12/2019 9:14 AM	File folder	
App 4 Name	None	A	Desktop	📮 SIGNU	IP	15/02/2023 6:02 PM	File folder	
			Documents	💷 ExtExp	port	07/12/2019 9:09 AM	Application	54 k
			Email attac	hm 🧟 iediag	ıcmd	07/12/2019 9:09 AM	Application	503 k
App 5 Name	None	A	Music	🦲 ieinsta	al	07/12/2019 9:09 AM	Application	493 k
			iviusic	🦲 ielow	util	07/12/2019 9:09 AM	Application	220 k
			Pictures	🤶 iexplo	re	12/03/2022 12:15 AM	Application	815 k
			Nandini Chav	vda				
		_	S This PC					
			🧮 Libraries					
			🐓 Network	~ <				>
				File name: iexplore			Executable Files (*.e	exe) ~
							Open	Cancel

FISECARE App Settings	3			×		
Example: App name: B	edside Navigator	App File Path:	C:\Program Files\Bedside Navigator\Bedside Navigat	or.exe		
Enter the name	Enter the name of the app		Select the absolute path of the app starting with C:\			
App 1 Name	Internet	App 1 File Path	C:\Program Files\Internet Explorer\iexplore.exe	Select File		
App 2 Name	None	App 2 File Path	None	Select File		
App 3 Name	None	App 3 File Path	None	Select File		
App 4 Name	None	App 4 File Path	None	Select File		
App 5 Name	None	App 5 File Path	None	Select File		
			Save	Cancel		

Then, save the file paths.

Step 4: A window will pop up to show what the camera is capturing. **Adjust the position of the camera until your face can be detected** well by MotionInput. A message will show to say when your face is detected saying 'Face Detected'.

Please wait a few seconds and stay in position until the calibration window shows.

Step 5: Below is a screenshot of how the **calibration window** looks like for face angle calibration.

The user should follow the dots by pointing to it with their nose.



Otherwise, for hand calibration, the user should follow the dots using their hand.

Information for calibration:

- Please wait for calibration to start. A black screen with instructions will show.
- Dots will show to guide your face or hand.
- Hold your hand or face at the position whilst a dot flashes.
- Don't move too far you should feel comfortable and still be able to see the instructions on the screen.
- Use the 'ESC' key on your keyboard to exit calibration whilst it is running.

Step 6: Once the dot reaches the bottom-right corner, calibration will end. The window will close and MotionInput will run as normal in a new window.

Setting Up (requires Keyboard and Mouse)

Step 1: Open the FISECARE folder and then open the Settings folder.

Step 2: In the **sender_email.json** file, enter SMTP host name, your care home email, the email's password, and the port number.

Save the file before closing. Please ensure you follow the format below. Alternatively, you could use the settings provided below:



Step 3: Go to <u>https://myaccount.google.com/</u> and log in to the email address you entered in step 2. Click on 'Security' on the left side and scroll down to 'Less secure app access'. **Please ensure this option is on.** This is how it should look like:



Step 4: In the **carer_email.txt** file, enter the carer's email and **save** the file. **Only one email is allowed**. Alternatively, you could use the provided email:

Carer_email -	Notepad		_		\times	
File Edit Format	t View Help					
carersoffise@gmail.com						
					~	
<					2	
Ln 1, Col 1	100%	Windows (CRLF)	UTF	-8 with B	OM	

The password for the email above is "Password_123".

Step 5: In the **family_email.txt** file, enter emails of residents' family members **(one email per line)** and **save** the file.

🧾 family_email - Not	tepad		_		\times	
File Edit Format Vie	ew Help					
familymember1@gmail.com						
familymember2@gmail.com						
familymember3@g	gmail.c	om				
< >						~
In 3 Col 24	100%	Windows (CRLE)	LITE-8		M	
LII 5, COI 24	10070	Windows (CILL)	011-0		VI.	

Step 6: In the **websites.json** file, you could change the URL for each more or leave it as it is. **Please ensure that you follow the format below and do not delete any existing modes.**

🧾 websites - Notepad × File Edit Format View Help { "Browser": "https://www.google.com", "Library": "https://manybooks.net/" "Call Cabs": "https://www.uber.com" "Grocery Delivery": "https://www.deliveroo.com", "Online Shopping": "https://www.amazon.com", "Watch Party": "https://www.netflix.com", "Classes": "https://www.youtube.com/channel/UCtFRv902AHq0Zjjynzrv-xg", "Drawing": "https://sketch.io/sketchpad/", "Games": "https://www.google.com", "Karaoke": "https://www.karafun.com", "Music": "https://open.spotify.com", "Chat": "https://web.whatsapp.com", "Video": "https://www.skype.com/en/free-conference-call/", "Calculator": "https://www.desmos.com/scientific", "Scratch": "https://scratch.mit.edu/", "Paint": "https://jspaint.app/#local:9c8333c2105d/", "WordLadder": "https://www.turtlediary.com/game/word-ladder.html", "Scrabble": "https://playscrabble.com/", "Study": "https://www.ixl.com/", "UNO": "https://buddyboardgames.com/uno", "MathsGames": "https://toytheater.com/category/math-games/", "HandyMaths": "https://mas281.github.io/hackathon/", "Sri Guru Granth Sahib": "https://www.srigranth.org/servlet/gurbani.gurbani?S=y", "Bhagavad Gita": "https://bhagavadgita.io/", "Tripitaka": "https://tripitaka.online/sutta/17", "Tanakh": "https://www.sefaria.org/texts", "Bible": "https://bible.com", "Quran": "https://quran.com", "Jaina Sutras": "https://www.sacred-texts.com/jai/", "ManyBooks": "https://freekidsbooks.org/" Ln 1, Col 1 100% Unix (LF) UTF-8

Login to websites (Keyboard and mouse required)

Step 1: Launch FISECARE. You may need to install <u>.NET Desktop</u> <u>Runtime 6.0</u>

Image: Image of the second	Manage Share View Application Tools	FISECARE				- 0	× ^ Ø
Pin to Quick Copy I access Clip	Paste Shortcut	py Delete Rename	New item •	Properties • Open • Open	Select all Select none Invert selection Select		
← → • ↑]	> This PC > Documents > FISE	CARE	5 v	, Search FISECA	RE		
🖈 Quick access	Name	~	Date modified 4/1/2022 10:45 PI	Type M File folde	Siz	e	^
OneDrive	Pictures		4/1/2022 10:45 PI	M File folde	r.		
This PC	runtimes		4/1/2022 10:45 PI	M File folde	f		
Network	FISECARE.deps		4/1/2022 10:45 PI 4/1/2022 10:45 PI	M File folde M JSON File	, 	33 KB	
	FISECARE.dll		4/1/2022 10:45 PI	M Applicatio	on extension	1,150 KB	
	fisecare		4/1/2022 10:45 PI	M Applicatio	n	163 KB	
	FISECARE.pdb		4/1/2022 10:45 PI	M Program	Debug Data	25 KB	
	FISECARE.runtimecon	fig	4/1/2022 10:45 PI	M JSON File		1 KB	
	S IronPython.dll		1/17/2022 2:24 PI	M Applicatio	on extension	1,773 KB	
	IronPython.Modules.c	111	1/17/2022 2:24 PI	M Applicatio	on extension	784 KB	
	IronPython.SQLite.dll		1/17/2022 3:59 PI	M Applicatio	on extension	623 KB	
	IronPython.Wpf.dll		1/17/2022 3:59 PI	M Applicatio	on extension	11 KB	
	Microsoft.Dynamic.dll	t,	12/2/2021 2:36 A	M Applicatio	on extension	811 KB	
28 items 1 item	selected 163 KB	e i e - 19 - 19 - 19 - 19		a a a a a a a a a a a a a a a a a a a	25. 1		

Step 2: Click on the **Services** button.



Step 3: Log in is required in the following pages:

- Under Concierge:
 - \circ Call Cabs
 - Grocery Delivery
 - \circ Online Shopping
- Under Creativity
 - o Karaoke
 - o Music
 - Watch Party

Services





Step 4: Select the page by clicking on the button.

Step 5: Click on the Login button.



Step 6: Login to the website.

	Music	
$\leftarrow \rightarrow$ https://accounts.sp	otify.com/login?continue=https%3A%2	Q
	Spotify [®]	
	To continue, log in to Spotify. CONTINUE WITH FACEBOOK CONTINUE WITH APPLE CONTINUE WITH GOOGLE	
	Email address or username Password Password	Pinay - Tan
← Creativity	MotionInput OFF E	TA FF

Step 7: Navigate back to the Home Page.

Step 8: Click on the 'For Kids' Button.

IBM in	tel. Augu	UCL FISECARE	Your Apps	Exit
		Wed, 22 March 2023 9:54 PM		
C	Exercises	Services	For Kids	
Created By IB	M			



Step 9: Click on the 'Scratch Coding' Button.

Step 10: If you don't already have an account for

https://scratch.mit.edu/ then skip to Step 12. Otherwise, click on the 'Sign in' button.

Scratch			
$\leftarrow \rightarrow$ https://scratch.mit.edu/		Ð	Q
Create Explore Ideas About Q Search Join Scratch	Sign in		Î
🛠 Scratch is the world's largest free coding community for kids. Your support makes a difference.			*
Create stories, games, and animations Share with others around the world Start Creating Join Watch Video	6		
← Kids Page MotionInput OFF	KIT/ OF	A F	

Step 11: Enter your username and password and sign in.



Step 12: Click on the 'Join Scratch' Button.

Scratch			
\leftarrow \rightarrow https://scratch.mit.edu/			@ Q
Create Explore Ideas About Q Search	Join Scratch	Sign in	Í
Scratch is the world's largest free coding community for kids. Your support makes a difference.	Donate		۲
Create stories, games, and animations Share with others around the world * Start Creating Join	Watch Vide		
← Kids Page OFF	Ţ	KITA OFF	

Step 13: Enter your details and follow the direction on the webpage to complete your registration.

Scratch											
4	→ https://scratch.mit.edu/join					Ð	Q				
(ALL III)											
			Join Scratch Create projects, share ideas, make friends. It's free! Create a username Username Create a password Password Type password again Not password again								
Next											
	←	- Kids Page	MotionInput OFF	₽	KITA OFF						

Navigating FISECARE with MotionInput

These steps show how you may navigate FISECARE to access a library of kids' books.

Step 1: Click on the **'For Kids'** Button.



Step 2: Click on the 'Online Books' button.

For Kids							
Scratch Coding	Calculator	Paint Word La	adder Online Books				
Scrabble	e UNO	Handy Maths	Study				
← Home							

Step 3: To **zoom in/out** of webpages, click the buttons on the top right corner. On the top left corner are the **Go Back /Forward** buttons to navigate you forward or back webpages. To **scroll up/down**, use the **scrollbar** on the right side of the screen to move up/down.

If using MotionInput's **speech commands** then say **'page up'** or **'page down'** to scroll up/down.

If using MotionInput's hand mode then you can use your index finger and middle fingers put together to scroll up/down.

If using MotionInput's face mode, open your mouth and face downward or upwards to scroll up/down



Step 4: Click on **'Read Online'** to read the book on the webpage OR click **'Download PDF'** to open the book in a new window as a PDF file.



Step 5: Click on the **'Speech To Text'** button to turn on the **speech-to-text** mode. Then, click on the search box and speak to search for the item. Here, we search for a book. Turn the mode off by clicking on the 'Speech To Text' button again.

