





User Manual

Discover how to harness the power of eye gaze with MI3.4 Eye Gaze. Here's a refined guide to get you started:

1. Launch the application by following the steps outlined in our deployment manual.
2. A Multiple Function Controller (MFC) dialog will prompt you to select your preferred settings.



3. A start-up window will appear. Before moving forward, make sure you're seated comfortably.
4. Next, a heat map will display to check if you're positioned at an optimal distance from the camera. A green indicator means you're perfectly positioned, while a red indicator suggests you need to adjust your distance.



5. Once you're correctly situated, the calibration process begins, displaying nine points on your screen sequentially.

Focus on each point as it appears.



6. After completing the calibration, the window will close, and you're all set to navigate using just your eyes!

Enjoy the seamless integration of MI3.4 Eye Gaze into your digital experience.
